
SAMPLE DATA WORKSHEET FOR MERCURY POROSIMETRY

_____ Basic presentation sample file

_____ Advanced presentation sample file

	Port 1	Port 2	Port 3	Port 4
Sample file name or identifier <i>(optional)</i>				
Penetrometer number <i>(optional - etched on penetrometer)</i>				
Sample mass(g) . Enter this value on the <i>Sample Description</i> tab.				
Sample + penetrometer mass . Load the penetrometer with sample and weigh.				
Penetrometer mass . Subtract the <i>sample mass</i> from the <i>Sample + penetrometer mass</i> and enter the difference. <ul style="list-style-type: none">• If using Basic presentation, enter this value on the <i>Sample Information</i> window.• If using Advanced presentation, enter this value on the <i>Penetrometer Properties</i> window.				
High pressure port number <i>(optional)</i> . Enter the number of the high pressure port where the sample is loaded.				
Assembly mass . The mass of the <i>sample + penetrometer + mercury</i> . (Weigh the penetrometer after the low pressure analysis.) Enter this value on the <i>Start High Pressure Analysis</i> window.				