

# Penetrometer Properties Worksheet for Mercury Porosimetry

File Name: \_\_\_\_\_

<b>Penetrometer Number</b> (etched on penetrometer)	
<b>Mass.</b> Weigh the loaded penetrometer. Subtract the mass of the sample. Enter the difference.	
<b>Volume</b>	
<b>Constant</b>	
<b>Stem Volume</b>	
<b>Max. head pressure</b>	
<b>Correction method</b> <input type="checkbox"/> Blank <input type="checkbox"/> Formula <input type="checkbox"/> None	If Blank, give file name: